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Type de document : imprimé divers

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1. Description

The Maltese cross tube is used to demonstrate the linear movements of cathode rays in a vacuum.

The apparatus consists of a conical glass bulb with two electrodes, the disk-shaped cathode at the conical end of the tube and the pointed anode fused in an appendage at the bottom.

A hinged cross made of aluminium sheet is used as the obstacle for the electrons flowing linearly through the bulb. The cross casts a sharply outlined shadow on the fluorescing face of the bulb. The pressure prevailing in the tube is approximately 0.01 mm. Hg.

An induction coil (521 11 or 521 02) is used as the source of voltage. To operate it, the tube is mounted with its appendage on the wooden stand supplied with the apparatus.

2. Operation

The tube is connected to the correctly poled induction coil; the cathode to the plate and the anode to the tip (for correct poling of the induction coil see Directions for Use 521 02 or 521 11).

When the tube is connected up correctly, a dark cross on a light base immediately appears on the face of the bulb.

Notes

1. The five-figure numbers quoted in brackets refer to the catalogue numbers of the respective apparatus.
2. The data DK . . . refer to the experiment descriptions published as a collection of leaflets in the "LEYBOLD PHYSIKALISCHE HANDBLÄTTER".
3. The data and illustrations are not binding in every detail for the design of the apparatuses. It is our sole aim always to adapt our manufacturing programme to the most recent knowledge gained in all scientific and technical fields.

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